**DE****Gebrauchsinformation**

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MEDICOPLAST INTERNATIONAL GmbH
Heusweilerstrasse 100
D-66557 Illingen

CE 1254

Instructions for Use
Seldinger Pleura Drainage Set with kink-resistant catheter

General information
REF 850 Seldinger Pleura Drainage Set with kink-resistant catheter (12 FR), latex-free, suitable for MRI

Indications
The Seldinger pleura catheter is used for draining serous, haemorrhagic serous or haemorrhagic pleural effusion; to relieve a pneumothorax.

Contraindications
The Seldinger pleura catheter is not suitable for draining a hemothorax with coagulated parts. When applying the set to patients with coagulation disorders or patients on anticoagulant medication, a thorough benefit-risk-analysis must be carried out before the application because of the potential risk of bleeding.

Instructions for Use
Before the puncture is carried out, make sure that there is actually an effusion or a pneumothorax underneath the selected puncture point and in the chosen puncture direction. In order to determine the puncture point and puncture direction, sonography is the method of choice with an effusion. In a pneumothorax, the puncture point and direction may be determined using a CT or thorax X-ray.

After a local anaesthetic has been applied to the puncture area, a small stab incision is made at the puncture point. This incision should be small in order to ensure that the skin tightly fits the catheter on all sides later on, so that the risk of bleeding and infection is minimised. If the skin section is extended with the scalpel when the drainage is in situ, there is the danger of damaging the drainage.

Now, with the tip applied, use the Tuohy needle (enclosed in the set) to puncture under continuous aspiration. The Tuohy needle enters the tissue in a way whereby the tissue is displaced. This means that the anatomical structure that is being penetrated can be felt while doing so. The puncture should be made at the top edge of a rib in order to avoid damaging the intercostal vessels that run alongside the bottom edge of the ribs. The scale on the needle allows you to reliably determine the current puncture depth which can be compared to the puncture depth that has been determined sonographically or radiologically.

In order to further increase the safety of the puncture, the puncture can be made under direct sonographic visualisation of the tip of the needle. After fluid or air has been aspirated, the needle should only be pushed forward for a very small distance and aligned in the direction where the wire is to run. Aspire once again to check that the tip of the needle continuously to be in the right position. Then, advance the Seldinger wire via the needle. After the applied tip has been removed and before the Seldinger wire is inserted into the cannula, make sure that no air can enter the thorax of the spontaneously breathing patient. This can be done, for example, by closing the cannula's opening with the thumb of the fixed hand. When the wire hits a spring-like resistance, at the very latest, the wire must not be advanced any further.

This process must be done while checking the heart rhythm with a monitoring ECG because the advancing wire can trigger arrhythmias. As soon as the wire is in situ, remove the Tuohy

needle, but make sure that the wire remains in its position. Then, the access route is dilated to the necessary diameter with the dilator. In order to do this, hold the dilator a short distance above skin level and, while slightly rotating it, insert it to just the depth that is necessary for a dilation down to the pleura. While doing this, check at regular intervals that the wire runs freely inside the dilator and has not been kinked while the dilator has been advanced.

After the dilation has been completed, remove the dilator and instead insert the drainage catheter into the thorax through the wire to a depth that the catheter's tip is positioned in the required area. Now, remove the Seldinger wire together with the inner stylet of the catheter. Screw the enclosed 3-way tap (in its closed position) onto the catheter. For this changeover from the Seldinger wire to the 3-way tap, you must again make sure that no air is able to enter the thorax via the drainage catheter.

The fluid can be drained into the drainage bag (enclosed in the set) by aspirating it with the enclosed tip. Alternatively, the drainage can be connected to a commercially available drainage system using the enclosed adapter. When a pneumothorax is relieved, the drainage must be connected to a drainage system. Now if required, the position of the drainage tip in the effusion or pneumothorax should be optimised by pulling the catheter back. This is because the catheter, if positioned too deeply, may become attached to the pleural gap. The depth marking on the catheter prevents the catheter from being pulled out too far. Once the perfect position has been found, you can use the enclosed plaster to fix the catheter or alternatively sew it into position. In order to check the correct position of the catheter and the success of the treatment, an overview X-ray of the thorax a.p. should be taken.

The catheter may be removed if the daily drainage quantity is below 300 ml with an effusion. It should be removed if the daily quantity is below 200 ml or the incision point shows a clear reddening. When relieving a pneumothorax the catheter should be removed when the fistula volume is reliable at zero millilitre. When removing the drainage, ensure that the tension is even, do not pull it in a jerky manner.

The maximum time of the catheter in situ is 28 days.

Storage
Store the catheter in a dry location and protected from direct and indirect light, moisture and heat sources.

Delivery
Each set is delivered packed individually in a sterile condition in a sleeve. The sterility of the set and its components is only guaranteed while the packaging unit remains unopened and undamaged, but at the most until the sterility date indicated on the label has expired.

Disposable products
Reliable cleaning and sterilisation methods as well as instructions for reconditioning according to the original specifications are not available for these products. The reuse of disposable products may result in their reduced performance, the loss of their correct functioning or in their exposure to pathogens. Validated verification procedures are not available.

Instructions for Use
Seldinger Pleura Drainage Set
with kink-resistant catheter (12 FR)

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TRACOE medical GmbH