EN Product Description - Hans Rudolph 7600 Series V2 Mask - Oro-Nasal CPAP/BiLevel Reusable - With Vents and AAV

Intended Use: The 7600 V2 Mask is a reusable, multi-patient, multi-use, Oro-Nasal CPAP/BiLevel mask which incorporates a passive, continuous flow exhaust port. It is intended for use on adults (>30 kg) in homes, hospitals and other clinical settings with CPAP/BiLevel machines for obstructive sleep apnea treatment, and for use with other ventilators that use this mask exhaust port configuration providing minimum 3 cm H2O pressure at mask.

Caution: Federal law restricts this device to be used by or on physician's order.

- Warnings:**
- 1. Do not use clean, nonsterile. If disinfected or sterile use is required, follow the procedures in this document prior to use of the mask.
 - 2. To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than a mild detergent.
 - 3. Patients with facial hair may experience mask leakage.
 - 4. Do not block mask vent holes or patient circuit flow.
 - 5. Do not use with ventilators that require double-limb patient circuits.
 - 6. Do not use with CPAP/BiLevel devices which include separate patient circuit exhaust port.
 - 7. When CPAP/BiLevel device is turned off, or is not operating properly, flow through mask vent holes may be inadequate to clear all exhaled gas and some rebreathing may occur which could lead to suffocation.
 - 8. Discomfort or irritation or allergic reaction may develop due to mask use.
 - 9. Patient failure to assemble and remove mask after vomiting could result in vomitus aspiration.
 - 10. If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
 - 11. Oxygen flow must be turned off when CPAP/BiLevel system is not operating.
 - 12. Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.

Contraindications:

- 1. Open wounds prone to infection.
- 2. Hemodynamic or cardiorespiratory instability.
- 3. Unconsciousness or patients unable to remove mask.
- 4. Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
- 5. Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
- 6. Excessive reflux, GI bleed or other secretions.
- 7. Impaired cough reflex, nasal hernia, or inability to swallow or clear secretions.
- 8. Discomfort with medication with respiratory or facial trauma.
- 9. Need for ventilation or ventilator support > 12 hours per day.
- 10. Recent facial, esophageal, or gastric surgery.
- 11. Patients under medication with respiratory or facial trauma.
- 12. Patients requiring immediate intubation.

Risks and Complications:

- 1. Biosensitivity to device materials (irritation, sensitivity, or other reaction).
- 2. Irritation/abrasion/corrosion/pain/deformity/itching/redness/eyes, mucosa, dental).
- 3. Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals).
- 4. Infection (skin, eyes, mucosa, respiratory).
- 5. Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma).
- 6. Excessive/damaging sound level.
- 7. Low circuit resistance especially during upper respiratory infections).

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High Level Disinfection (Mask and Swivel port components see Headgear)

- 1. Thoroughly clean mask prior to disinfection or sterilization.
- 2. For disinfection use Code™ liquid glutaraldehyde or hot water pasteurization. Liquid Chemical Disinfection: Soak the device components in the liquid chemical solution per the manufacturer's instructions. Remove the components and soak in 1500 ml of sterile water for a minimum of 1 minute. Repeat step 2 and then dry the device with a clean lint free cloth (preferably sterile).
- 3. For hot water pasteurization immerse all device components surfaces in hot water temperatures 71.1°C – 76.6°C for 30 minutes. Dry the device with a clean lint free cloth (preferably sterile).
- 3. Visually inspect all components for cleanliness, function and defects. If there are any residues, stains or organic debris that are not the previous steps. Dispose of and replace all defective parts.

Reassemble:

- 1. Snap AAV into swivel port end, ensuring both elbow locking tabs are fully engaged with slots of the AAV.
- 2. Press rubber cap plug onto sampling port of mask adapter.
- 3. Install the mask adapter assembly into face mask groove opening. The flange of the mask adapter mates with the groove of the mask opening.
- 4. Mount headgear to face piece with the snap-in quick-release clips.

Functional Check:

- 1. Turn on oxygen supply freely; swivel port assembly is completely engaged in mask face.
- 2. AAV is installed completely in port with locking tabs fully engaged and the diaphragm flexes freely with no obstruction and functions as described in the safety features section. Do not use if it is non-functional.

Storage:

- 1. Components should be completely dry prior to storage in a sealed poly bag.
- 2. Label poly bag with disinfection/clean status, mask description, date and initials.

Steam Sterilization (Face Mask assembly less Headgear)

- 1. Température: 132.2 +/-3/1 °C
- 2. Sterilization time: 4 minutes
- 3. Dry time: 10 minutes
- 4. Precaution: No use solutions based on chlorine or a base of peroxy or alcohol to clean the mask. If the mask is touched or wrapped in CSR.

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