	ESPAÑOL (es)	DEUTSCH (DE)	FRANÇAIS (FR)
A	ADAPTADOR MASCARILLA	MASKENADAPTER	ADAPTEUR DE MASQUE
B	PIEZA FRONTAL	GESICHTSTEL	PIECE FACIALE
C	TOMA MUESTREO	PROBENNAHMEANSCHLUSS	ORIFICE D'ÉCHANTILLONNAGE
D	VALVULA AAV	AAV	VALVULE VAA
E	LENGÜETAS TRABANTES	ENRÄSTZUNGEN	LANGUETTES DE VERROUILLAGE
F	TOMA GIRATORIA	DREHGELENKSCHLUSS	RACCORDS ARTICULÉS
G	ACODADA	GEBÖGENER TYP	TYPE COUDÉE
H	CONEXIONES COFIA	KOPFHÄUTER-BEFESTIGUNG	FENTES POUR ATTACHES DU HARNAIS

	ESPAÑOL (es)	DEUTSCH (DE)	FRANÇAIS (FR)
A	LENGÜETAS TRABANTES	ENRÄSTZUNGEN	LANGUETTES DE VERROUILLAGE
B	PRENSION HASTA QUE QUEDEN TRABADAS EN SU LUGAR	HINEINDRÜCKEN BIS ES EINSCHNAPPT	ENFONCEZ JUSQU'À CE LA POSITION S'OT VERROUILLÉE
C	VALVULA AAV	AAV	VALVULE VAA

	ESPAÑOL (es)	DEUTSCH (DE)	FRANÇAIS (FR)
A	TOMA AAV	AAV-ANSCHLUSS	ORIFICE VAA
B	DIAPHRAGMA	PLATTENVENTIL	DIAPHRAGME
C	DIAPHRAGMA	PLATTENVENTIL	DIAPHRAGME
D	ABIERTA	OFFEN	OUVERT
E	CERRADA	GESCHLOSSEN	FERMÉ

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EN Product Description - Hans Rudolph 7700 Series V2 Mask - Oro-Nasal BiLevel Reusable - With AAV and NO Vents

Intended Use: The 7700 V2 Mask is a reusable, multi-patient, multi-use, Oro-Nasal BiLevel mask that does NOT incorporate a passive, continuous flow exhaust port. It is intended for use on adults (>30 kg) in homes, hospitals and other clinical settings with CPAP/BiLevel machines for obstructive sleep apnea treatment, and for use with other similar ventilators that incorporate an exhaust port into the patient circuit instead of the mask and provide minimum 3 cm H₂O pressure at mask.

Caution: Federal law restricts this device to sale by or on physician's order.

- Warnings:**
- Spigged chain, nonsterile. If disinfected or sterile use is required, follow the procedures in this document prior to use of the mask.
 - To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than mild detergent and water.
 - Patients with facial hair may experience mask leakage.
 - Do not block mask or patient circuit flow.
 - Do not use with ventilators that require double-limb patient circuits.
 - Not for use with CPAP/BiLevel devices which require a mask with built-in vent holes.
 - When BiLevel device is turned off, or is not operating properly, flow through mask and vented patient circuit may be inadequate to clear all exhaled gas and some rebreathing may occur which could lead to suffocation.
 - Discontinue mask use if irritation or allergic reaction develops due to mask.
 - Patient failure to arouse and remove mask after vomiting could result in vomiting aspiration.
 - If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
 - Oxygen flow must be turned off when CPAP/BiLevel system is not operating.
 - Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.

- Contraindications:**
- Open wounds prone to infection.
 - Hemodynamic or cardiorespiratory instability.
 - Unc consciousness or patients unable to remove mask.
 - Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
 - Facial (or nasopharyngeal) deformity or other inability to fit mask & seal properly.
 - Excessive reflux, GI blood or other secretions.
 - Impaired cough reflex, nasal hernia, or inability to swallow or clear secretions.
 - Barotrauma, upper airway obstruction or facial trauma.
 - Need for ventilation or ventilatory support > 12 hours per day.
 - Recent facial, esophageal, or gastric surgery.
 - Patients under medication with drug that may cause vomiting.
 - Patients requiring immediate intubation.

- Risks and Complications:**
- Bio-sensitivity to device materials (irritation, sensitivity, or other reaction).
 - Inhalation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental).
 - Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals).
 - Infection (skin, eyes, mucosa, respiratory).
 - Excessive/rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma).
 - Excessive/damaging sound level.
 - Low secretion clearance especially during upper respiratory infections).
 - Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence).
 - Adverse interactions with attached non-HRI accessories.
- Safety Information:**
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- Mask Components (fig 1):**
- Mask face piece has plastic bracing with slot openings for attachment of headgear strap mounting quick-release clips.
 - Swivel port assembly comprises sampling port, anti-asphyxia valve (AAV) and detachable 22 mm (ISO 5356-1) swivel port.
 - AAV (fig 2) snaps into swivel port assembly elbow using "locking tabs" for engagement.
- Application & Fitting (fig 4):**
- Determine patient mask size with optional Mask Sizing Gauge (fig 3).
 - Place over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
 - Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
 - Adjust headgear Velcro™ strap tension by pulling back slightly on both the top and bottom straps (4 & 5). Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
 - Connect mask to patient circuit by following ventilator operating instructions.
 - Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.
- Safety Features Verification:**
- AAV (fig 5) allows room air breathing for patient if ventilator device stops. AAV diaphragm CLOSES (>3 cm H2O) the room air port in mask during CPAP/BiLevel ventilation. It OPENS (<3 cm H2O) mask to room air breathing when ventilator device is off. Mask should not be used if AAV is missing, damaged or not functioning properly.

- Mask Removal with Quick-Release Headgear (fig 6). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

- Disassembly for Cleaning**
- Separate headgear from mask.
 - Separate swivel port assembly from mask. AAV from swivel port (squeeze two locking tabs & pull AAV out) and rubber cap plug from mask adaptor sampling port.

- Cleaning Mask face piece, swivel port components & headgear.**
- Mask & swivel port assembly:**
- Soak components for 5 minutes in warm water with mild detergent (neutral pH).
 - Hand wash components with sponge or soft brush; then rinse with warm tap water.
 - Place components in warm water bath and agitate for two minutes.
 - Rinse in clean tap water for at least 1 minute.
 - Allow components to air dry or dry with clean, lint free cloth.

- Headgear:**
- Submerge headgear in warm soapy tap water; gently rub all areas.
 - Rinse in clean tap water for 1 minute or until all soap is removed; air dry.
 - Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions can damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.

- High Level Disinfection (Mask and Swivel port components less Headgear)**
- Thoroughly clean mask prior to disinfection or sterilization.
 - For disinfection use Cide™ liquid glutaraldehyde or hot water pasteurization.


- Liquid Chemical Disinfection**
- Soak the device components in the liquid chemical solution per the manufacturer's instructions. Remove the components and soak in 1500 ml of sterile water for a minimum of 1 minute. Repeat step 2 and then dry with a clean lint free cloth (preferably sterile).

- Hot Water Pasteurization**
- Immerse all device components surfaces in hot water temperatures 71.1°C - 76.6°C for 30 minutes. Dry the device with a clean lint free cloth (preferably sterile).

- Visually inspect all components for cleanliness, function and defects.** If there are any residues, stains or organic debris then repeat the previous steps. Dispose of and replace all defective parts.

- Reassembly:**
- Snip AAV into swivel port elbow, ensuring both elbow locking tabs are fully engaged with slots of the AAV.
 - Press rubber cap plug onto sampling port of mask adaptor.
 - Install swivel port assembly into face mask groove opening. The flange of the mask adaptor mates with the groove of the mask opening.
 - Mount headgear to face piece with the snap in quick-release clips.

- Functional Check:**
- Swivel port joints swivel freely, swivel port assembly is completely engaged in mask face piece.
 - AAV is installed completely in port with locking tabs fully engaged and its diaphragm flexing freely with no obstruction

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