
FIGURE 1		
ESPAÑOL (ES)	DEUTSCH (DE)	FRANÇAIS (FR)
A ADAPTADOR MASCARILLA	MASKENADAPTER	ADAPTATEUR DE MASQUE
B PIEZA FRONTAL	GESICHTSTEIL	PIÈCE FACIALE
C TOMA MUESTREO	PROBENABNAHMENSCHLUSS	ORIFICE D'ÉCHANTILLONNAGE
D TOMA GIRATORIA	DREHGELEKENSCHLUSS	RACCORD ARTICULÉ
E RECTA	GERADER TYP	MODÈLE DROIT
F TOMA GIRATORIA	DREHGELEKENSCHLUSS	RACCORD ARTICULÉ
G ACODADA	GEBOGENER TYP	TYPE COUDÉ
H CONEXIÓN COFIA	KOPFHAUBEN-BEFESTIGUNG	FENTES POUR ATTACHES DU HARNAIS

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Intended Use:
The 6500 V2 Mask is a disposable single-patient use, Oro-Nasal NIV mask intended to provide a patient interface for applications of noninvasive ventilation. The mask is intended for use on adults (>30 kg) in hospitals and other clinical settings as an accessory to ventilators which have adequate alarms and safety systems for ventilator failure and which are intended to administer positive pressure ventilation for treatment of respiratory failure or respiratory insufficiency.

Caution:
Federal law restricts this device to sale by or on physician's order.

Warnings:

- Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions.
- To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than mild detergent and water.
- Inspect and clean before use, replace if not cleanable.
- Patients with facial hair may experience mask leakage.
- Do not block mask or patient circuit flow.
- This mask does not include an Anti-Asphyxia Valve (AAV) that would allow patient to breathe if ventilator failed. Use only with ventilation device that has anti-suffocation mechanism and adequate alarms and safety systems for ventilator failure.
- Do not use with single-limb constant-flow ventilation or CPAP devices that require an exhalation valve, vent holes or anti-suffocation valve built into the mask.
- Discontinue use if patient has allergic reaction develops due to mask.
- Patient failure to breathe and remove mask after vomiting could result in vomitus aspiration.
- If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
- Oxygen flow must be turned off when the ventilation system is not operating.
- Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.

Contraindications:

- Open wounds prone to infection.
- Hemodynamic or cardiorespiratory instability.
- Unconsciousness or patients unable to remove mask.
- Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
- Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
- Excessive reflux, GI blood or other secretions.
- Impaired cough reflex, hiatal hernia, or inability to swallow or clear secretions.
- Barotrauma, upper airway obstruction or facial trauma.
- Recent facial, esophageal, or gastric surgery.
- Patients under medication with drug that may cause vomiting.
- Patients requiring immediate intubation.

Risks and Complications:

- Biosensitivity to device materials (irritation, sensitivity, or other reaction).
- Irritation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental).
- Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals).
- Infection (skin, eyes, mucosa, respiratory).
- Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma).
- Excessive/damaging sound level.
- Low secretion clearance (especially during upper respiratory infections).
- Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence).
- Adverse interactions with attached non-ARI accessories.

Safety Information:
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Mask Components (Fig 1):

- Mask face piece has plastic braces with slot openings for attachment of headgear strap mounting quick-release clips.
- Swivel port assembly (elbow or Straight Type) comprises sampling port, 22 mm (ISO 5356-1) swivel port.

Application & Fitting (Fig 3):

- Determine patient mask size with optional Mask Sizing Gauge (Fig 2).
- Place mask over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
- Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
- Adjust headgear Velcro™ strap tension by pulling back slightly on the both top and bottom straps (4 & 5). Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
- Connect mask to patient circuit by following ventilator operating instructions.
- Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.

Safety Features Verification:

- Mask Removal with Quick-Release Headgear (Fig 4). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

Disassembly for Cleaning

- Separate headgear from mask.
- Separate swivel port assembly from mask and rubber cap plug from mask adapter sampling port.

Cleaning Mask face piece, swivel port components & headgear.

- Soak components for 5 minutes in warm water with mild detergent (neutral pH).
- Hand wash components with sponge or soft brush; then rinse with warm tap water.
- Place components in warm water bath and agitate for two minutes.
- Rinse in clean tap water for at least 1 minute.
- Allow components to air dry or dry with clean, lint free cloth.

Headgear:

- Submerge headgear in warm soapy tap water; gently rub all areas.
- Rinse in clean tap water for 1 minute or until all soap is removed; air dry.
- Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions can damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.

Reassemble:

- Press rubber cap plug onto sampling port of mask adapter.
- Install swivel port assembly into face mask grooved opening.
- Mount headgear to face piece with the snap in quick-release clips.

Functional Check:

- Swivel port joints swivel freely; swivel port assembly is completely engaged in mask face piece. Do not use if it is non-functional.

Technical Specifications

- Operating pressure: 0 – 40 cm H₂O
- Mask Resistance to flow: 0.3 cm H₂O at 50 L/min, 0.7 cm H₂O at 100 L/min
- Unintentional mask leakage (leak at swivel port connections)

Pressure (cm H ₂ O)	10	20	30
Leak (L/min)	3	4	5

- Deadspace volume (ml)

L(169)	M(151)	S(125)	XS(114)	P(104)
10	20	30	40	50

Mask & Headgear Service Life:
1. Mask and headgear assembly may be cleaned and reused on single-patient basis 7 days max. Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions.

Operational Environment, Mask Assembly & Headgear:

- Operational Temperature: 5 – 40 °C, Humidity: 0 – 95% RH
- Storage Temperature: -40 – 70 °C, Humidity: 0 – 95% RH

Recommendations for Disposal:

- Treat as conventional solid waste in accordance with local and federal regulations.

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