	ESPAÑOL (es)	DEUTSCH (de)	FRANÇAIS (fr)
A	ADAPTADOR MASCARILLA	MASKENADAPTER	ADAPTATEUR DE MASQUE
B	PIEZA FRONTAL	GESICHTSTEIL	PIÈCE FACIALE
C	TOMA MUESTRO	PROBENABNAHMANSCHLUSS	ORIFICE D'ÉCHANTILLONNAGE
D	ORIFICIOS VENTILACION	LUFTBOHRE	ORIFICES D'ÉRATION
E	VALVULA AAV	AAV	VALVULE VAA
F	LENGUETAS TRABANTES	EINRASTZUNGEN	LANGUETTES DE VERROUILLAGE
G	TOMA GIRATORIA	DREHGELENKANSCHLUSS	RACCORDES ARTICULES
H	ACCODADA	GEBOGENER TYP	TYPE COUDÉ
I	CONEXIONES COFIA	KOPFHAUBEN-BEFESTIGUNG	FENTES POUR ATTACHES DU HARNAIS

	ESPAÑOL (es)	DEUTSCH (de)	FRANÇAIS (fr)
A	LENGUETAS TRABANTES	EINRASTZUNGEN	LANGUETTES DE VERROUILLAGE
B	PRESIONE HASTA QUE QUEDEN TRABADAS EN SU LUGAR	HINEINDRÜCKEN BIS ES EINSCHNAPPT	ENFONCEZ JUSQU'À CE LA POSITION SOIT VERROUILLÉE
C	VALVULA AAV	AAV	VALVULE VAA

	ESPAÑOL (es)	DEUTSCH (de)	FRANÇAIS (fr)
A	TOMA AAV	AAV-ANSCHLUSS	ORIFICE VAA
B	DIAPHRAGMA	PLATTENVENTIL	DIAPHRAGME
C	DIAPHRAGMA	PLATTENVENTIL	DIAPHRAGME
D	ABIERTA	OFFEN	OUVERT
E	CERRADA	GESCHLOSSEN	FERMÉ

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HANS RUDOLPH, INC.
6860 Series QUEST™ Mask
Oro-Nasal CPAP
BiLevel
with Vents & AAV
Single Patient Use
Home Care Setting Only
Instructions for Use

EC REP E C Rep Ltd.
9 Fitzwilliam Square East
Dublin 2, D02 R744
Ireland
Tel: +353 1 2 544 944
Email: info@ecrep.ie

PRINTED IN U.S.A.
691284 2021-06-11 Rev. E

HANS RUDOLPH, INC.
8325 Cole Parkway
Shawnee, Kansas 66227 USA
Tel: (913) 422-7788
Toll Free: (800) 456-6695
Fax: (913) 422-3337
Web Site: www.rudolphk.com
Email: hr@rudolphk.com

Medical Device

MD

US Pat No's D573,709 & 8,136,523 and
EC Req'd Design No. 00078056-001

© 2797 © 2021 Hans Rudolph, inc.

EN Product Description – Hans Rudolph 6860 Series Quest™ Mask – Oro-Nasal CPAP/BiLevel – Single Patient Use – Home Care – With Vents and AAV

Intended Use:
The 6860 Quest™ Mask is a single-patient use, Oro-Nasal CPAP/BiLevel mask which incorporates a passive, continuous flow exhaust port. It is intended for use on adults (>30 kg) in the home care setting only with CPAP/BiLevel machines for obstructive sleep apnea treatment, and for use with other similar ventilators that use this mask exhaust port configuration providing minimum 3 cm H₂O pressure at mask.

Warnings:
Federal law restricts this device to sale by or on physician's order.
Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions. To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than mild detergent and water.
Inspect and clean before use, replace if not cleanable.
Patients with facial hair may experience mask leakage.
Do not block mask vent holes or patient circuit flow.
Do not use with ventilators that require double-limb patient circuits.
Not for use with CPAP/BiLevel devices which include separate patient circuit exhaust port.
When CPAP/BiLevel device is turned off, or is not operating properly, flow through mask vent holes may be inadequate to clear all exhaled gas and some re-breathing may occur which could lead to suffocation.
Discontinue mask use if irritation or allergic reaction develops due to mask.
Patient failure to arouse and remove mask after vomiting could result in vomitus aspiration.
If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
Oxygen flow must be turned off when CPAP/BiLevel system is not operating.
Combination of this mask with devices or other than for which mask is intended can alter mask fit and performance becoming unsafe or life-threatening.

- Contraindications:**
- Open wounds prone to infection.
 - Hemodynamic or cardiorespiratory instability.
 - Unconsciousness or patients unable to remove mask.
 - Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
 - Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
 - Excessive reflux, GI blood or other secretions.
 - Impaired cough reflex, hiatal hernia, or inability to swallow or clear secretions.
 - Barotrauma, upper airway obstruction or facial trauma.
 - Need for ventilation or ventilatory support > 12 hours per day.
 - Recent facial, esophageal, or gastric surgery.
 - Patients under medication with drug that may cause vomiting.
 - Patients requiring immediate intubation.

Risks and Complications:

- Biosensitivity to device materials (irritation, sensitivity, or other reaction).
- Irritation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental).
- Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals).
- Excessive skin, eyes, mucosa, respiratory).
- Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma).
- Excessive/damaging sound level.
- Low secretion clearance (especially during upper respiratory infections).
- Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence).
- Adverse interactions with attached non-HRI accessories.

Safety Information:
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Mask Components (Fig 1):

- Mask face piece has plastic braces with slot openings for attachment of headgear strap mounting quick-release clips.
- Swivel Port assembly comprises sampling port, micro vent holes, anti-asphyxia valve (AAV) & detachable 22 mm (ISO 5356-1) swivel port.
- AAV (Fig 2) snaps to the swivel port assembly elbow using "locking tabs" for engagement.

Application & Fitting (Fig 4):

- Determine patient mask size with optional Mask Sizing Gauge (Fig 3).
- Snep mask over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
- Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
- Adjust headgear Velcro™ strap tension by pulling back slightly on the top and bottom straps (4 & 5). Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
- Connect mask to patient circuit by following ventilator operating instructions.
- Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.

Safety Features Verification:

- AAV (fig 5) allows room air breathing for patient if ventilator device stops. AAV diaphragm CLOSÉS (>3 cm H₂O) the room air port in mask during CPAP/BiLevel ventilation. It OPENS (<3 cm H₂O) mask to room air breathing when ventilator device is off. Mask should not be used if AAV is missing, damaged or not functioning properly.
- Mask Removal with Quick-Release Headgear (Fig 6). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

Disassembly for Cleaning:

- Separate headgear from mask assembly.
- Soak mask assembly in warm water with mild detergent for 2 minutes. The warm water with detergent assists in disconnecting the swivel port assembly from the mask.
- Separate the swivel port assembly from mask by squeezing the face piece close to the mask adapter while pulling mask adapter out at 90 degrees.
- Separate the AAV from the swivel port by squeezing the two locking tabs and pull the AAV out. Remove the rubber cap plug from mask adapter sampling port.

Cleaning Mask face piece, swivel port components & headgear

Mask & swivel port assembly:

- Soak components for 5 minutes in warm water with mild detergent (neutral pH).
- Hand wash components with sponge or soft brush; then rinse with warm tap water.
- Place components in warm water bath and agitate for two minutes.
- Rinse in clean tap water for at least 1 minute.
- Allow components to air dry or dry with clean, lint free cloth.

Headgear:

- Submerge headgear in warm soapy tap water; gently rub all areas.
- Rinse in clean tap water for 1 minute or until all soap is removed; air dry.
- Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions may damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.

Reassembly:

- Snap AAV into swivel port elbow, ensuring both elbow locking tabs are fully engaged with slots of AAV.
- Press rubber cap plug onto sampling port of mask adapter.
- Install swivel port assembly into face mask grooved opening.
- Mount headgear to face piece with the snap in quick-release clips.

Functional Check:

- Swivel port joints swivel freely; swivel port assembly is completely engaged in mask face piece.
- AAV is installed completely in port with locking tabs fully engaged and its diaphragm flexing freely with no obstruction and functions as described in the safety features section. Do not use it if it is non-functional.

Technical Specifications

Pressure Flow Chart for Mask Exhaust Vent

Operating pressure: 4 – 40 cm H₂O
Mask Resistance to flow: 0.5 cm H₂O at 50 L/min, 1.0 cm H₂O at 100 L/min
AAV Resistance to flow (open port): 3.5 cm H₂O at 50 L/min
AAV opening/closing pressure: 1.5 – 3 cm H₂O
Unintentional mask leakage (leak at swivel port connections):

Pressure (cm H ₂ O)	10	20	30
Leak (L/min)	5	8	10

Deadspace volume (ml):

L154	M136	S110	XS(99)	P(89)
------	------	------	--------	-------

Sound pressure level: <25 dB(A) measured 15 cm from mask

Mask & Headgear Safety Life:
Mask and headgear assembly may be cleaned and reused on single-patient basis for 6 months use under normal conditions. Do not attempt to sterilize or disinfect mask or headgear. Sealing surfaces and headgear will lose shape under these processes causing excessive leaks. Follow cleaning instructions.

Operational Environment, Mask Assembly & Headgear:

- Operational: Temperature: 5 – 40 °C, Humidity: 0 – 95% RH
- Storage: Temperature: 40 – 70 °C, Humidity: 0 – 95% RH

Recommendations for Disposal:

- Treat as conventional solid waste in accordance with local and federal regulations.

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Installation and Functional Check:

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