FIGURE 1								
	ESPAÑOL (ES)	DEUTSCH (DE)	FRANCAIS (FR)	ITALIANO (IT)	NEDERLANDS (NT)	SVENSKA (SV)	PORTUGUES (PT)	ARABIC (AR)
A	Toma hembra de 22 mm	22 mm Schlauchanschluss	Orifice femelle 22 mm	Porta femmina da 22 mm	22 mm Vrouwelijke Poort	22 mm honkoppeling	Porta fêmea de 22 mm	منفذ أنثى قيس 22 مم
B	Presillas de desenganche rápido de la cofia	Kopfhäube Quick-Release	Attaches rapides du harnais	Rilasci rapidi della cuffia	Snelkoppelings Hoofdband	Headset snab-koppoling	Liberações rápidas do headgear	إعتاق سريع لعدة الرأس
C	Conexiones de la cofia	Kopfhäubenbefestigungsstellen	Fente pour attache du harnais	Accessori della cuffia	Toevoegings Hoofdband	Headset fästen	Acessórios do headgear	مُلحقات عدة الرأس
D	Toma hembra de muestra Luer	Luer Probenabnahmschluss	Orifice d'échantillonnage Luer femelle	Porta di campionamento per luer femmina	Vrouwelijke Luer Koppelingpoort	Provtagningsport, Luer-honkoppeling	Porta de amostragem luer fêmea	منفذ عينات Luer أنثى
E	Toma de sellado endoscópico	Endoskopisch abdichtender Anschluss	Orifice hermétique pour endoscopie	Porta di sigillatura endoscopica	Endoscopische Afdichtingpoort	Tätslutande port för endoskop	Porta de vedação endoscópica	منفذ الختم المنظار
F	Toma de muestra de espigas de manguera	Schlauchtilen-Probenabnahmschluss	Orifice d'échantillonnage à raccord cannelé	Porta di campionamento per portagomma	Koppeling-spool Slang met Haak	Provtagningsport för rafflad slangkoppeling	Porta de amostragem da espiga da mangueira	منفذ عينات قرنة الخرطوم

HANS RUDOLPH, INC.

DISPOSABLE 7302 SERIES

EAGLE 2™ FACE MASK

INSTRUCTIONS FOR USE

E C Rep Ltd.
5 Fitzwilliam Square East
Dublin 2, D02 R744
Ireland
Tel: +353 1 2 544 944
Email: info@ecrep.ie

PRINTED IN U.S.A.
691902 2021-06-11 Rev. D

HANS RUDOLPH, INC.
8325 Cole Parkway
Shawnee, Kansas 66227 USA
Tel: (913) 422-7788
Toll Free: (800) 456-6695
Fax: (913) 422-3337
Web Site: www.rudolphkc.com
Email: hri@rudolphkc.com

2797 © 2021 Hans Rudolph, Inc.

US Pat D 761, 416 and EC
RCD: 002661884

EN HANS RUDOLPH, INC. 7302 SERIES EAGLE 2™ FACE MASK INSTRUCTIONS FOR USE

Intended Use:
The 7302 Eagle 2 Mask is a disposable single-patient use, Oro-Nasal mask intended to administer oxygen and other breathing gases into the upper airway for short-term or emergency applications of noninvasive respiratory support or ventilation in hospital, transport vehicle, institutional or other clinical settings on patients that weigh > 30 kilograms. The Mask has a supplemental, membrane-sealed port for use during insertion of a bronchoscope or other endoscopic device into the patient's airway while still administering oxygen and other breathing gases.

- Caution:**
Federal law restricts this device for use ONLY by medical professionals.
- Warnings:**
1. Disposable Single-Patient Use only.
 2. Not for Home Care use.
 3. Do not attempt to sterilize or disinfect mask or headgear. They will lose shape under these processes causing excessive leaks. When necessary, follow cleaning instructions.
 4. Patients with facial hair may experience mask leakage.
 5. This mask does not include anti-asphyxia valve (AAV) that allows patient to breathe if ventilator fails. Use a ventilator with built-in anti-suffocation mechanism and adequate alarms and safety systems for ventilator failure.
 6. Use of this mask for CPAP or BiLevel applications requires an external vent valve or port in patient circuit or ventilator.
 7. Discontinue mask or headgear use if irritation or allergic reaction develops.
 8. Patient failure to arouse and remove mask after vomiting could result in vomitus aspiration.
 9. If supplemental oxygen is used, inhaled percent oxygen could vary depending upon pressure settings, patient breathing pattern, mask size and leak rate.
 10. Oxygen flow must be turned off when the ventilation system is not operating.
 11. Combination of this mask with devices other than for which mask is intended can alter its safety and performance, becoming unsafe or life-threatening.

- Contraindications:**
1. Open wounds in mask area prone to infection;
 2. Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask;

ES HANS RUDOLPH, INC. MASCARILLA FACIAL EAGLE 2™ SERIE 7302 INSTRUCCIONES DE USO

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3. Facial or nasopharyngeal deformity or other inability to fit mask and seal properly;
4. Excessive reflux, GI blood or other secretions;
5. Impaired cough reflex or inability to swallow or clear secretions; (unless medically induced during endoscopy procedure);
6. Barotrauma, upper airway obstruction or facial trauma;
7. Recent facial, esophageal, or gastric surgery;
8. Patients under medication with drug that may cause vomiting; and
9. Patients requiring immediate intubation.

- Risks and Complications:**
1. Biosensitivity to device materials (irritation, sensitivity, or other reaction);
 2. Irritation/abrasion/ulceration/pain/deformity/drying (face, head, neck, eyes, mucosa, dental);
 3. Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals);
 4. Infection (skin, eyes, mucosa, respiratory);
 5. Excessive re-breathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma);
 6. Excessive/damaging sound level;
 7. Low secretion clearance (especially during upper respiratory infections);
 8. Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence); and
 9. Adverse interactions with attached non-HRI accessories.

Safety Information:
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

- Mask Components (fig. 1):**
1. Mask face piece has slit openings for attachment of headgear straps.
 2. Headgear straps have quick-release clips for removing mask from patient's face
 3. Main mask port is 22 mm female connection per (ISO 5356-1)
 4. Hose barb and female Luer sampling ports
 5. Endoscopic sealing port

- Application and Fitting (fig. 3):**
1. Determine patient mask size with sizing gauge instruction printed on the mask bag (fig 2). With the patient's mouth closed, position the gauge as shown in Fig 2 for best mask size selection.
 2. With the headgear attached to the mask (as packaged), place the headgear over the patient's head while positioning mask on the patient's face (1).

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3. Fit mask chin cup with patient's chin, then move top of mask sealing area onto the bridge of the nose (2). Slightly rock mask left-to-right and top-to-bottom to assure good fit and position on the face.
4. While holding mask on the patient's face, adjust each of the headgear straps by pulling each strap's plastic clip back and away engaging the hook & loop until the headgear is tight and adjusted equally (3). Be sure to position top of headgear open crown portion on top of head for stabilization.
5. Repeat the mask, strap and headgear adjustments until mask is positioned properly and secure on the patient's face.

Safety Features Verification:
Mask removal with quick-release headgear: Disengage lower headgear strap's hook & loop by pulling on the strap's plastic clips, then pulling mask or headgear over the patient's head. The straps do not have to disengage the mask for complete removal from the patient.

Cleaning of Soiled Mask:
Use warm soapy water with a good rinse and dry, if needed. Remove tethered caps when cleaning and drying. Do not scrub or distort endoscopy membrane seal during cleaning process. Reset caps and seal ports after drying (as required).

Functional Check:
Check tethered cap setting and seal for all three supplementary ports

- Technical Specifications:**
1. Operating pressure: 0 – 40 cm H2O.
 2. Deadspace volume (ml): L (145), M (110), and S (80).
 3. Endoscopic instrument diameters accepted: 3.6 to 6.3 mm

Mask and Headgear Service Life:
Under the normal range of operational and environmental conditions, as defined below, the mask face piece and headgear are expected to stay in service for a minimum of 24 hours of operation. These components are disposable after use and are incapable of being disinfected, sterilized without compromising their performance.

- Environment (Mask and Headgear):**
1. Operational: -18 to 50 °C; 0 to 100% RH
 2. Extreme Ambient Transport Conditions (Temporary): -40 to 70 °C; 0 – 95% RH

Recommendations for Disposal:
Treat as conventional solid waste in accordance with local and federal regulations.

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